# Leeds-Ruff Patch™

**EN** Instructions for Use

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







## Description

The Leeds-Kuff Patch is a permanent implantable device for reinforcement of rotator cuff tears. It is constructed from knitted fabric with integral reinforcements across its width and around its perimeter to provide improved stability and high suture retention strength. The Leeds-Kuff Patch does not require cutting as it is available in three sizes that are pre-shaped to match typical rotator cuff tears. The patches taper over their length and are approximately 20 mm x 20 mm, 30 mm x 30 mm, and 40 mm x 35 mm.

## Material Specifications

The Leeds-Kuff Patch is made from polyethylene terephthalate (polyester), which meets applicable specifications established by ISO 10993-1 for implantable medical devices.

## Intended Use

The Leeds-Kuff Patch is intended for reinforcement of soft tissues, where weakness exists, that are repaired by suture or suture anchors during rotator cuff repair surgery. The Leeds-Kuff Patch is not intended to replace normal body structure or provide the full mechanical strength to support the rotator cuff. The Leeds-Kuff patch reinforces soft tissue and provides a scaffold that is incorporated in the patient's own tissue.

## Indications

The Leeds-Kuff Patch is indicated for patients requiring reinforcement of the rotator cuff where either the tear cannot be completely repaired using normal methods and / or the quality of the soft tissue is poor.

## Contraindications

- Known hypersensitivity to implant materials. If the patient is suspected of having any foreign body sensitivity, appropriate tests should be made prior to implantation.
- Infections or any structural or pathological condition of the bone or soft tissue that would be expected to impair healing or secure fixation.
- Patients unable or unwilling to restrict activities to prescribed levels or follow a rehabilitation programme during the healing period.
- Any repair that requires rolling, folding, or layering, and which may create a space impermeable to fluid, cells, and blood vessels. Such uses may result in excessive inflammation, drainage, extrusion or infection.
- The Leeds-Kuff Patch should not be used to provide full mechanical strength to the rotator cuff as it is intended to reinforce tears that are repaired with sutures and / or bone anchors.
- The Leeds-Kuff Patch integrates well with soft tissue so may not be suitable where eventual removal of the device is anticipated.

## Warnings

- The device is provided sterile and remains so unless the packaging is damaged or opened. Store in standard conditions. DO NOT RESTERILIZE.
- The device is for single use only as it would be damaged by reprocessing. Do not use after the expiration date. Discard any open, unused product.
- The surgeon must be thoroughly familiar with these instructions and the recommended surgical procedure before using the device.
- The general principles of patient selection and sound surgical judgement apply to the reconstruction procedure.
- The Leeds-Kuff Patch must not be cut to size.

## Precautions

- Inspect the device, packaging and labelling prior to use and do not use if damaged. Contact the Neoligaments Sales department if the package has been opened or altered.
- Take care to select the most suitable size of Leeds-Kuff Patch for the surgery to be undertaken, based on patient anatomy, size of tear, or extent of poor quality tissue.

- If using an arthroscopic technique use a minimum 8 mm cannula.
- All sutures attached for load transfer must be placed through the patch, inside and adjacent to the reinforced perimeter (Figure 1). Load bearing sutures must not be placed through or outside the reinforcement, or towards the centre of the patch, as they may break or pull through the device. Supplementary sutures used to ensure apposition of the patch to the underlying cuff may be placed within the centre of the patch.
- Simple stitches are preferred for attaching to tissue as they are easier to perform in restricted joint space and provide high suture retention strength when encircling the perimeter reinforcement. Alternatively mattress stitches may be used.
- When using suture anchors or per-osseous sutures take care to leave an adequate bone thickness between the tunnel and the bone surface or between adjacent tunnels. Take the quality of the bone into account.
- Check that any knots tied in the sutures are securely locked.
- Avoid damage when handling the Leeds-Kuff Patch. Avoid crushing or crimping when using surgical instruments such as forceps or needle holders.
- It is recommended that wherever possible patients are placed on prophylactic antibiotics prior to surgery, to minimize the risk of latent infections developing at the implant site.
- Use aseptic technique throughout the procedure.
- Patients should be warned not to exceed appropriate activity levels or to overload the repair before complete healing has occurred.

### Potential Adverse Effects

- Wound dehiscence, infected wounds, minimal acute inflammatory tissue reaction, and transitory local irritation. If any of these conditions occur, the device may need to be removed at the surgeon's discretion.
- The sutures used to secure the Leeds-Kuff Patch may cause elongation of the device, or the sutures may break through the soft tissue, if prescribed activity levels during rehabilitation are exceeded.
- Tissue erosion.

### Instructions for Use: Open Approach

1. With the patient in the beach chair position use a lateral deltoid splitting approach to achieve access to a posteriorly subluxed infraspinatus. Dissect the anterior-lateral deltoid subperiosteally in the line of its fibres to a minimum of 5 cm from

the lateral edge of the acromion. The lower limit of this dissection is the axillary nerve. A stay suture can be applied across the deltoid to prevent further splitting.

2. Gently apply a self-retaining retractor and perform a subacromial bursectomy if necessary to improve vision. A subacromial decompression may also be performed if joint space is limited, taking care to avoid damaging the suprascapular nerve.
3. Distract the shoulder and rotate to access the anterior and posterior aspects of the cuff repair. Place #2 high strength stay sutures in the cuff, leaving the needles attached. If access is limited #1 suture may be used.
4. Repair the defect in the rotator cuff as far as possible, tying the stay sutures for a side to side repair where achievable (Figure 2). The tension should be moderate at most.
5. Select the most suitable size of Leeds-Kuff Patch for the repair based on patient anatomy, size of tear or extent of poor quality tissue.
6. Create simple stitches with #2 high strength suture to attach the edges of the patch over the superior edges of the rotator cuff tendon so that it sits flush with the tendon. Ensure load bearing sutures are placed through the patch, inside the inner border of the perimeter reinforcement. Avoid over- or under-tensioning.
7. Start with a vertical stitch at the apex then place stitches at the medial edge anteriorly and posteriorly (Figure 3). Place further stitches on the anterior and posterior edges of the patch. If the size of tear is large or the quality of the tendon is weak, more than one stitch on each side may be required to distribute the load.
8. Use the stay sutures placed previously to provide supplementary attachment to ensure apposition of the patch to the underlying cuff. As these sutures are non-load bearing they may be placed within the centre of the patch.
9. Pull the repair laterally and use #5 suture with an appropriate integral needle to anchor the patch and tendon to the tuberosity (Figure 4). Up to three per-osseous sutures may be needed, depending on the size of the patch. Take care to leave sufficient bone both between adjacent sutures and the bone surface to resist expected forces. Take into account the quality of the bone.
10. Alternatively, suture anchors with high strength suture can be used in place of the per-osseous sutures.
11. Apply appropriate tension to ensure the length of the reconstruction is physiological. Avoid over- or under-tensioning to ensure adequate joint function is achieved and premature failure of the repair does not occur.

- 
- Carefully repair the deltoid proximally with mattress stitches. Repair the edges of the incision in the deltoid. An appropriate suture with a stitch for fat can be used if applicable and a subcuticular suture to skin.

### **Instructions for Use: Arthroscopic Technique**

- With the patient in the beach chair or lateral decubitus position apply gentle traction to distract the subacromial space to facilitate viewing.
- Make a standard posterior viewing port referencing from the postero-lateral corner of the acromion. Enter the subacromial space and under direct vision use a needle to identify the lateral port to allow optimal direct vision of the cuff tear.
- Repair the defect in the rotator cuff as far as possible with moderate tension at most. Leave the sutures from the suture anchors used in the repair for later fixation of the anterior and posterior edges of the patch (Figure 5). For the small patch use these for medial corner fixation (step 6).
- Select the most suitable size of Leeds-Kuff Patch for the repair based on patient anatomy, size of tear or extent of poor quality tissue. For arthroscopic insertion and fixation the patch is positioned with the apex oriented laterally.
- Use the Neoligaments Leeds-Punch (supplied separately) to pierce holes in the patch in the preferred locations for sutures. Ensure load bearing sutures are placed through the patch, inside and adjacent to the reinforced perimeter (Figure 1). Generally the small patch will only need openings at the medial corners and the lateral apex.
- Prior to insertion of the patch place two #2 high strength polyester sutures through the rotator cuff such that once tied to the patch one will form the antero-medial and one the postero-medial corner suture fixation.
- Retrieve one limb of each of the passed sutures through the lateral cannula and tie to the respective anterior or posterior corner of the medial border of the patch, taking care not to twist the sutures. The free limbs of the passed sutures are taken out through accessory antero-medial and postero-medial ports respectively (Figure 6).
- Fold or pleat the patch along its length and alternately tension the free ends of the sutures to pull the patch into the subacromial space via the lateral cannula (Figure 7).
- Manipulate the patch into the required position and tie the free ends of the sutures (the post suture) to the corresponding suture previously tied to the corners of the patch. Avoid over- or under-tensioning (Figure 8).
- Place a suture anchor laterally. Pass sutures from it through the lateral edge of the patch and tension.
- Apply appropriate tension to ensure the length of the reconstruction is physiological. Avoid over- or under-tensioning to ensure adequate joint function is achieved and premature failure of the repair does not occur.
- For medium and large patches use the sutures from the previously placed suture anchors to secure the anterior and posterior edges of the patch to the underlying cuff. Alternatively, separate stay sutures may be used. If the size of tear is large or the quality of the soft tissue is weak, more than one stitch on each side may be required to distribute the load.
- Add an additional medial suture if required.
- Close the lateral port with an appropriate interrupted skin suture.

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All rights reserved. © Neoligaments™ 2016.
Leeds-Kuff Patch, Neoligaments and Xiros are trademarks of Xiros.
Xiros Limited, registered in England No. 1664824.