neoligaments™

# JewelACL™

**EN** Instructions for Use

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


102-6003 JewelACL™, 7 mm ID x 710 mm

### Description

The JewelACL is a polyester prosthesis for anterior cruciate ligament (ACL) reconstruction. It is matched in tensile strength to the semitendinosus hamstring tendon, which is typically 1200 N, and can be used with or without an additional tissue graft in either partial or total tissue sparing ACL reconstruction procedures.

### Material Specifications

The JewelACL is made from polyethylene terephthalate (polyester). The polyester has been subjected to a proprietary gas plasma process that modifies its surface properties, making it hydrophilic, without significantly altering the physical characteristics of the bulk material. This increases the speed and extent of cell recruitment and adhesion and the remodelling of the tissue into a “neoligament”.

### Fixation

The JewelACL must only be secured using the approved fixation devices detailed below. When using these devices it is critical to follow the instructions for use supplied with them and to observe the conditions and necessary adjustments stated in the notes below.

### Femoral Fixation

EndoButton® CL Ultra (Smith & Nephew) OR  
TransFix® (Arthrex)

### Tibial Fixation with a Screw

Interference Screw, 7 mm x 30 mm (Medgal) 5.5 mm  
bone tunnel, OR  
RCI Screw®, 7 mm x 25 mm (Smith & Nephew)  
6.5 mm bone tunnel

### Tibial Fixation with a Screw and Tissue Graft

Size the bone tunnel to match the overall graft diameter and use a larger diameter screw. For example:

Graft diameter	6.5 mm
Tunnel diameter	6.5 mm
Screw diameter	8.0 mm

This applies to normal bone density with a typical graft diameter. At the discretion of the surgeon the screw diameter may be increased in cases encountering soft bone in order to obtain secure fixation.

### Intended Use

The JewelACL is intended to be used for reconstruction of the anterior cruciate ligament.

### Indications

The JewelACL is indicated for all patients requiring ACL reconstruction that are not excluded by the contraindications listed below.

### Contraindications

- Known hypersensitivity to implant materials. If the patient is suspected of having any foreign body sensitivity, appropriate tests should be made prior to implantation.
- Infections, or any structural or pathological condition of the bone or soft tissue that would be expected to impair healing or secure fixation.
- Patients unable or unwilling to restrict activities to prescribed levels or follow a rehabilitation programme during the healing period.
- The JewelACL may not be suitable for skeletally immature patients as it will not elongate and so must not bridge, disturb, or disrupt the growth plate.
- Patients for whom it is not possible to bend the knee to at least 90° as it will not be possible to reach the correct position for drilling the bone tunnels.

### Warnings

- The device is provided sterile and remains so unless the packaging is damaged or opened. Store in standard conditions between 5°C and 30°C. DO NOT RESTERILIZE.
- The device is for single use only as it would be damaged by reprocessing. Do not use after the expiration date. Discard any open, unused product.
- The surgeon must be thoroughly familiar with these instructions and the recommended surgical procedure before using the device.
- The general principles of patient selection and sound surgical judgement apply to the reconstruction procedure.
- There are no data to support the use of the JewelACL in pregnant and breast-feeding women.

### MRI Safety Information

The JewelACL is MR safe (i.e. an item that poses no known hazards in all MR environments).

### Precautions

- Inspect the device, packaging and labelling prior to use and do not use if damaged. Contact the Neoligaments Sales Department if the package has been opened or altered.
- Correct isometric placement is crucial to the successful outcome of the reconstruction. Failure to place the ligament bone tunnels correctly is

---

likely to lead to an unsatisfactory outcome. The ligament must be implanted taut while allowing full range of motion (ROM).

- Use of fixation devices other than those listed may damage the JewelACL and lead to potential failure of the device.
- Chamfer the edges of the bone tunnel exits where possible to minimize abrasion of the JewelACL.
- Ensure bone tunnel length is adequate when using an interference screw for tibial fixation. The screw must be centrally placed between the two strands of the JewelACL, which should then be tied over the back of the screw with a reef knot. Ensure the knot is locked before trimming any excess JewelACL.
- When using a tissue graft and an interference screw, ensure excess graft is trimmed from the tibial tunnel entrance prior to knotting the JewelACL.
- When trimming the JewelACL, cut straight across it at right angles to its length to minimize generation of loose fibre fragments. Carefully remove any loose fibres from the incision wound.
- Ensure that any knots remain well buried in tissue.
- Use aseptic technique throughout the procedure.
- Take care to avoid anatomical hazards, including soft tissue and nerve damage, when making incisions and drilling tunnels.
- Patients should be warned not to exceed appropriate activity levels before complete healing has occurred.
- Preoperative prophylactic antibiotics should be administered.

### Potential Adverse Effects

- Infections, both deep and superficial.
- Patient body sensitivity and adverse reaction to device materials.
- Laxity or reduced ROM are potential adverse consequences of misplacing the ligament.
- Laxity may also result from inadequate fixation of the ligament to bone.

### Disposal

No specific disposal requirements other than handling contaminated items as clinical waste.

### Instructions for Use

The JewelACL can be implanted using an anteromedial or transtibial surgical technique similar to that followed when reconstructing the ACL with hamstring grafts.

The following instructions describe an anteromedial approach incorporating a single hamstring tendon to produce a partial tissue sparing graft in which the JewelACL is fixed to the femur with an EndoButton CL Ultra (ECLU) and to the tibia with an interference screw.

1. Establish standard anterolateral (AL), anteromedial (AM) and accessory medial portals with the knee flexed at 90°.
2. Create a low medial accessory portal as low as possible but above the medial joint line while avoiding the anterior horn of the medial meniscus. The medial-lateral placement of this portal should be close to the medial edge of the patellar tendon to maximize the tunnel length for such a cortical suspension device.
3. Harvest an appropriate hamstring tendon and prepare according to standard procedures so that each end is whip-stitched with sutures to a length of approximately 35-40 mm.
4. Thread the ends of the sutures through the eyelet of a passing pin and pass it through the appropriate openings in the JewelACL, pulling the tendon inside. Do not pierce the JewelACL with the passing pin.
5. Where a shallow intercondylar notch is likely to cause impingement the knee should be flexed beyond 90° (typically 130°) and a femoral notchplasty performed through the AM portal.
6. Position an appropriate zero or small offset femoral aimer through the low medial accessory portal at the over-the-top position on the femur. The femoral tunnel exit in the intercondylar notch should be as far posterior and proximal as possible while remaining within the ACL footprint (Figure 1). Ensure that the tunnel has an adequate wall thickness to prevent breakthrough.
7. Drill a passing pin through the femoral aimer, into the femur and out through the anterolateral cortex. Avoid drilling into the peroneal nerve, or damaging the cartilage surface of the medial femoral condyle.
8. Remove the aimer, leaving the passing pin. Over-drill with a 4.5 mm cannulated bone tunnel drill until the cortex is breached, then measure the length of the tunnel.
9. Choose the appropriate size of ECLU and assemble with the graft following the technique described by the manufacturers to form a two-strand graft (Figure 2). Measure the graft diameter to determine the required bone tunnel size.
10. Drill the femoral socket to the appropriate depth with a drill matched to the graft size. Ensure the drill does not breach the lateral femoral cortex, or fixation with the ECLU cannot be performed. Where possible chamfer the edges of the distal tunnel at its intercondylar exit with an ACL tunnel rasp.
11. With the knee flexed at 90° use an appropriate tibial guide to identify the intra-articular tibial attachment (Figure 3). This should be located slightly medial and slightly anterior to the centre of attachment of the natural ACL. It should not interfere with the anterior attachment of the medial meniscus and should avoid damaging the articular cartilage. Avoid placing too anteriorly as this can lead to impingement of the ligament on the roof of

---

the notch at full extension. Leave the tibial footprint of the ACL intact.

12. Use a tibial guide to drill a guidewire into the tibia. Ensure an adequate tunnel length is produced to accommodate an interference screw.
13. Remove the tibial guide, leaving the guidewire. Over-drill with a cannulated bone tunnel drill matched to the overall graft size or make a tunnel 1 mm smaller than the graft and expand to the desired size using serial dilators. Where possible, chamfer the edges of the intercondylar exit with an ACL tunnel rasp.
14. Pull the assembly into the tunnels following the instructions for the ECLU. Ensure the button is securely seated on the femoral cortex.
15. With the knee in full extension, separate and tension the two strands of the graft. Introduce the screw centrally between the strands in the bone tunnel. Ensure the strands of the graft do not twist and end up on the same side of the screw.
16. Trim excess hamstring graft from the tibial tunnel entrance and tie the strands of the JewelACL over the back of the screw with a reef knot to provide additional fixation and prevent the graft slipping past the screw (Figure 4).
17. Where possible attach the ACL remnants to the graft using appropriate sutures.
18. Cycle the knee through a full range of motion while examining the graft arthroscopically to ensure that it has been placed isometrically and allows a full range of motion with no graft impingement.
19. Check security of the knot before trimming any excess strands of the JewelACL. Cut each strand with scissors at right angles to its length, to minimize the generation of loose fibres. Any loose fibres created when trimming to length must be carefully removed from the incision site.
20. Bury the knot in soft tissue.

### **Trademarks**

TransFix is a registered trademark of Arthrex, Inc.

EndoButton and RCI Screw are registered trademarks of Smith & Nephew, Inc.

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All rights reserved. © Neoligaments™ 2020.
trademarks of Xiros.
Xiros Limited, registered in England No. 1664824