# Poly-Tape

**EN** Instructions for Use

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## Description

Poly-Tapes are sterile non-absorbable implantable tapes made from polyester. Poly-Tapes are available for the Indications listed below as tubular or single thickness, in open or dense weave, in widths of 3-40 mm and 500-800 mm long (see Ordering Information). They are supplied for single use only.

## Material Specifications

The Poly-Tape is made from polyethylene terephthalate (polyester), which meets applicable specifications established by ISO 10993-1 for implantable medical devices.

## Intended Use

1. Poly-Tapes are single use devices intended to be used for tissue approximation, including use in reconstructing damaged or torn ligaments, tendons or other soft tissues according to the surgeon's own preferred technique, and at his/her discretion.
2. Poly-Tapes are also supplied for soft tissue (tendon and ligaments) fixation to bone with a Fastlok™ Fixation Device during orthopaedic reconstruction procedures.

## Indications

Poly-Tapes are indicated for patients requiring soft tissue approximation and reconstruction of ligaments, tendons and other soft tissues.

## Contraindications

1. Known hypersensitivity to implant materials. If the patient is suspected of having any foreign body sensitivity, appropriate tests should be made prior to implantation.
2. Infections or any structural or pathological condition of the bone or soft tissue that would be expected to impair healing or secure fixation.
3. Patients unable or unwilling to restrict activities to prescribed levels or follow a rehabilitation programme during the healing period.
4. The use of this device may not be suitable for patients who are skeletally immature. It will not elongate as the patient grows and so it must not bridge, disturb, or disrupt the growth plate. Bone quality should be assessed and surgery should not be performed on patients with insufficient or immature bone.
5. Poly-Tapes integrate well with the patient's tissue and consequently may not be suitable where eventual removal of the tape is anticipated.

## Warnings

1. The device is provided sterile and remains so unless the packaging is damaged or opened. Store in standard conditions. DO NOT RESTERILIZE.
2. The device is for single use only as it would be damaged by reprocessing. Do not use after the expiration date. Discard any open, unused product.
3. The general principles of patient selection and sound

surgical judgement apply to the reconstruction procedure.

4. Users should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing Poly-Tapes.
5. Take care to select a suitable Poly-Tape for the surgery to be undertaken.
6. Forming a knot in any tape or suture may reduce its overall strength by approximately 50%.
7. When used with a fixation device, the strength of the combination may be less than that of the tape used alone.
8. When used with a fixation device such as staples, using excessive force to impact the staple may cause damage to, or cut through, the Poly-Tape.

## Precautions

1. In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to the application of surgical instruments such as forceps or needle holders.
2. It is recommended that wherever possible patients are placed on prophylactic antibiotics prior to surgery, including oral surgery, to minimize the risk of latent infections developing at the device site.
3. If drilling bone tunnels to receive the Poly-Tape, take care to leave an adequate bone thickness between the tunnel and the surface or between 2 parallel tunnels, to provide a sufficient bony bridge to resist expected forces that will be exerted on the bone by the reconstruction. Take into account the quality of the bone.
4. Check that any knot is securely locked.
5. Be aware that the Poly-Tape might be damaged by sharp edges of bone tunnels and certain fixation devices. Wherever possible all bone edges should be rounded and fixation devices with rounded edges should be used.
6. To trim excess length after knotting, cut straight across the Poly-Tape at right angles to its edge to minimize generation of loose fibre fragments. After trimming it may be necessary to restrain the cut ends by stitching them back to the material.
7. Any loose fibres formed when trimming the Poly-Tape to length must be carefully removed from the incision wound.
8. Ensure that any knots remain well buried in tissue.
9. Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.
10. The patient should be warned not to exceed the prescribed activity levels or to overload the repair before complete healing has occurred.
11. Take care to avoid anatomical hazards including soft tissue and nerves. Ensure they are not trapped by the tape.

### Potential Adverse Effects


1. Possible adverse reactions associated with the use of this device include: wound dehiscence, infected wounds, minimal acute inflammatory tissue reaction, and transitory local irritation.
2. The tape could possibly break through the bone if bone tunnels are incorrectly positioned, the bone is of poor quality, or prescribed activity levels are exceeded.

### Directions for Use

Further information if required may be obtained from the Neoligaments Sales Department.

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All rights reserved. © Neoligaments™ 2018.
Fastlok, Neoligaments and Xiros are trademarks
Xiros Limited, Registered in England No. 1664824.