MPFL System Implant Set™

EN Instructions for Use

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

102-1060 MPFL System Implant Set, includes:

Poly-Tape, 10 mm x 500 mm
Fastlok™, 6 mm x 23 mm

Packaged with the following disposables:

Malleable probe with eye, stainless steel, 20 cm
Drill Bit, plain shank to fit Jacobs Chuck, 3.2 mm diameter

Description

The MPFL System comprises two permanent implantable devices that together enable reconstruction of the medial patellofemoral ligament in order to stabilize the patella against lateral translation. A 10 mm x 500 mm open weave Poly-Tape is looped through and around the patella and secured to the femoral attachment site with a 6 mm Fastlok staple, thus avoiding the need to harvest any autogenous tissue. The MPFL System is packaged with the required 3.2 mm drill bit and 20 cm malleable probe.

Material Specifications

The Poly-Tape is made from polyethylene terephthalate (polyester) and meets applicable specifications established by ISO 10993-1 for implantable medical devices. The Fastlok staple is made from titanium alloy (Ti-6Al-4V) and meets the ISO 5832-3 standard for metallic implants.

Intended Use

The MPFL System is intended to be used in the reconstruction of the medial patellofemoral ligament.

Indications

The MPFL System is indicated for patients with recurrent patellar subluxation or dislocation.

Contraindications

- Known hypersensitivity to implant materials. If the patient is suspected of having any foreign body sensitivity, appropriate tests should be made prior to implantation.
- Infections or any structural or pathological condition of the bone or soft tissue that would be expected to impair healing or secure fixation.

- Patients unable or unwilling to restrict activities to prescribed levels or follow a rehabilitation programme during the healing period.
- Skeletally immature patients. The MPFL System will not elongate with growth and must not bridge or disrupt the growth plate.
- Poly-Tapes integrate well with the patient's tissue and consequently may not be suitable where eventual removal of the tape is anticipated.

Warnings

- The contents of the MPFL System are provided sterile and remain so unless the packaging is damaged or opened. Store in standard conditions.
- DO NOT RESTERILIZE the Poly-Tape.
- The MPFL System is for single use only. The Poly-Tape would be damaged by reprocessing. The metallic items will be subject to wear during the technique which may affect performance on any subsequent use. Do not use after the expiration date. Discard any open, unused product.
- The surgeon must be thoroughly familiar with these instructions and the recommended surgical procedure before using the product.
- The general principles of patient selection and sound surgical judgement apply to the reconstruction procedure.

Precautions

- Inspect the devices, packaging and labelling prior to use and do not use if damaged. Contact the Neoligaments Sales Department if the package has been opened or altered.
- Care should be taken to avoid the saphenous nerve and other anatomical hazards during surgery.
- Care should be taken to diathermy branches of the descending geniculate and superior medial geniculate arteries at the femoral fixation site.
- Care should be taken to avoid tunnel blow-out when drilling the bone tunnel in the patella.
- Check that the repair is physiological and does not restrict range of motion (ROM).
- Care should be taken to avoid toggling the Fastlok when performing the ROM check.
- When trimming the Poly-Tape to length, cut straight across at right angles to minimize the generation of loose fibres and carefully remove any created from the wound.
- Avoid damage when handling the Poly-Tape. Avoid crushing or crimping when using surgical instruments such as forceps or needle holders.
- Aseptic technique should be followed throughout and wherever possible patients should be placed on prophylactic antibiotics prior to surgery to

minimize the risk of latent infections developing at the implant site.

- Patients should be warned not to exceed appropriate activity levels before complete healing has occurred.
- Caution is advised with patients presenting with either a Q angle greater than 25° or severe patella alta (patella height assessed by the Insall-Salvati index greater than 1.3).

Potential Adverse Effects

- Wound dehiscence, infected wounds, minimal acute inflammatory tissue reaction, and transitory local irritation.

Instructions for Use

The following provides instructions on how to use the MPFL System during a standard approach.

1. Make a 4-6 cm medial longitudinal incision two thirds of the distance between the medial epicondyle and the medial border of the patella to access the adductor tubercle and the lateral border of the patella (Figure 1).
2. Perform superficial dissection to the medial patellofemoral retinaculum. Take care to avoid damaging the infrapatellar branch of the saphenous nerve.
3. Mark the entry point for the patellar bone tunnel on the medial border of the patella within the centre of the MPFL footprint, approximately at the junction of the proximal and middle thirds of the patella (Figure 2). To avoid tunnel blow-out, take care not to place the tunnel too far anteriorly in the sagittal plane.
4. Use sharp dissection to raise a flap of the medial patellofemoral retinaculum from the point marked on the medial border of the patella to the point of origin of the ligament, just distal to the vastus medialis. The flap should be 10 mm wide by 6-8 cm long and should be left attached to the patella (Figure 2). If the synovial lining is breached it may be repaired using a 2/0 vicryl suture.
5. Locate the site of femoral fixation within the footprint of the origin of the MPFL, approximately equidistant from the medial epicondyle and the adductor tubercle (Figure 3). Take care to diathermy branches of the descending geniculate and superior medial geniculate arteries at the site.
6. To reduce the prominence of the Fastlok staple, use sharp dissection and then a Böhler bone rongeur to remove a 12 mm by 3 mm sliver of bone at the identified site, approximately 45° to the long axis of the femur (Figure 3).
7. With the knee in extension, use the supplied 3.2 mm drill bit to create a bone tunnel in the patella, passing from the identified point on the medial side (step 3) and emerging at the anterolateral border. If the exit point and orientation of the tunnel are not correct they may be repositioned by over-drilling with a 4.5 mm drill bit (available as an optional extra), which may also be used if a larger tunnel is preferred. Where possible, round the tunnel edges to prevent abrasion of the Poly-Tape.
8. Thread the Poly-Tape through the eyelet of the supplied probe and pass it through the bone tunnel from medial to lateral.
9. Pass Crile artery forceps across the anterior surface of the patella, deep to the anterior retinaculum, and pass the Poly-Tape back across the patella from lateral to medial.
10. Check the patella tracks well throughout the range of motion by holding the Poly-Tape to the attachment site while flexing and extending the knee. In order to avoid over-tightening, the Poly-Tape should be fixed with the knee positioned at the angle which creates the greatest tension. This is typically found in full flexion.
11. Using the technique supplied with the Fastlok, thread the free ends of the Poly-Tape onto the staple. Position the prongs at approximately 45° to the long axis of the femur, slightly forward towards the patella, avoiding entry into the knee joint.
12. Place Gillies forceps deep to the Poly-Tape during Fastlok insertion and apply a light tension to the free ends of the tape to ensure the Poly-Tape is not over-tightened. Confirm the correct position of the patella in the femoral groove.
13. Use the Fastlok Impactor instrument (sold separately) to impact the staple into the bone to a depth of approximately 1 cm while maintaining tension on the Poly-Tape. Perform a second check of the flexion and extension if required, to check the tension in the Poly-Tape. Release the Impactor from the staple before conducting this check to avoid toggling the staple. Check the tension in the Poly-Tape afterwards to ensure the tape did not slip past the Fastlok during the check.
14. If satisfactory, impact the staple as far as possible to ensure it is not palpable. Move the knee through a full range of flexion to ensure the Poly-Tape is not too tight.
15. Use scissors to trim any excess Poly-Tape to length, leaving a short tail from each end to place over the Fastlok. Cut the tape at right angles to its length to minimize the generation of loose fibres, and remove any that are created.
16. Replace the medial soft tissue flap and repair it along the superior and inferior edges to the surrounding retinaculum, suturing to the Poly-Tape wherever possible (Figure 4).
17. Suture through the capsule, Poly-Tape and retinaculum flap to anchor the cut ends of the tape. Take care to bury the ends of the tape and the staple in soft tissue (Figure 4).

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All rights reserved. © Neoligaments™ 2017.
are trademarks of Xiros.
Xiros Limited, registered in England No. 1664824.