neoligaments™

CC-Hook™

- EN** Instructions for Use
- ES** Instrucciones de uso
- DE** Hinweise zur Anwendung
- EL** Οδηγίες χρήσης
- FR** Consignes d'utilisation
- IT** Istruzioni per l'uso
- PT** Instruções de uso
- NL** Gebruiksaanwijzing
- SV** Användningsinstruktioner
- PL** Instrukcje użytkowania
- MS** Arahan penggunaan
- FI** Käyttöohjeet
- SL** Navodila za uporabo

---

CE 2797

LAB 229 4.00  
01 - 2020

## Ordering Information (supplied sterile)

202-1411 CC-Hook™ (left)

202-1413 CC-Hook™ (right)

### Patient Information

The following information is provided for use by clinicians, however as the learned intermediary between the company and the patient, the clinician must convey the aspects they consider relevant to the individual patient. The patient must be informed of the potential adverse effects (risks/complications) contained in this insert (see **POTENTIAL ADVERSE EFFECTS**).

### Description

The CC-Hook is a single-use handheld instrument with a curved end that is capable of passing under and around the coracoid process while causing minimal damage to the surrounding soft tissues. It is cannulated and supplied with a flexible wire that is passed through the length of the instrument so that it protrudes from the far end of the device to capture a synthetic ligament or suture which can then be pulled back under the coracoid process.

The CC-Hook is available in two variants, one for the left coracoid and one for the right. It is intended to be used by surgeons trained in acromioclavicular joint repair.

### Material Specifications

The CC-Hook comprises Stainless Steel, Polybutylene terephthalate (PBT) and Nitinol.

### Intended Use

The CC-Hook is intended to facilitate passing of synthetic ligaments or sutures under the coracoid process in acromioclavicular joint reconstruction.

### Indications For Use

The CC-Hook is indicated for use in patients undergoing acromioclavicular joint reconstruction with synthetic ligaments or sutures that are to be passed underneath the coracoid process.

### Contraindications

- The instrument must not be used for any procedure other than the intended use.

### Warnings

#### Pre-operative

- The device is single-use only as it is not possible to clean effectively. Do not use after the expiration date. Discard

any open, unused product.

- The surgeon must be thoroughly familiar with these instructions and the recommended surgical procedure before using the device.
- The general principles of patient selection and sound surgical judgement apply.

#### Intra-operative

- The CC-Hook is supplied with a protective cap. This must be removed prior to use.
- Caution should be exercised when handling devices with sharp points or cutting edges.

#### Precautions

##### Packaging

- The device is provided sterile (by gamma irradiation) and remains so unless the packaging is damaged or opened. Store in standard conditions. DO NOT RESTERILIZE.
- Inspect the device, packaging and labelling prior to use and do not use if damaged. Contact Xiros if the package has been opened or altered.

##### Handling and storage

- To avoid damaging the instrument do not impact it, subject it to excessive blunt force or use it to manipulate bone.
- Devices should be treated with care and improper use or handling may lead to damage and/or possible malfunction.
- Devices should be checked to ensure that they are in working order prior to surgery.
- All devices should be inspected prior to use to ensure that there is no unacceptable deterioration such as corrosion, discoloration etc. Non-working or damaged devices should not be used and should be returned to Xiros.
- No specific storage conditions are required other than good warehousing practice.

##### Pre-operative

- In morbidly obese patients it may be difficult to gain adequate access to the coracoid process.
- In unusually muscular patients it may be difficult to pass the instrument around the coracoid process.
- In unusually small patients it may be difficult to pass the instrument around the coracoid process without risking injury to the brachial plexus.

##### Intra-operative

- The instrument must not be used by untrained or unqualified personnel.
- Care must be taken to avoid nerve and circulatory system damage when passing the instrument under the coracoid process. The device should be passed in a medial to lateral direction.
- Use aseptic technique throughout the procedure.

---

## Potential Adverse Effects

Below is a list of the potential adverse effects (e.g. complications) associated with the use of the device including, 1) risks associated with any surgical procedure; 2) risks associated with coracoid passers, including the CC-Hook. Additional surgery may be required to correct some of these events.

- 1) Pertinent risks associated with any surgical procedure include: Infection and device breakage leaving material in situ.
- 2) Risks associated with coracoid passers, including the CC-Hook, include: nerve damage; numbness; paralysis; and surrounding tissue injury.

## Surgical Technique

1. Make a 4-5 cm incision centred over the coracoid ligament, approximately 3.5 cm medial to the acromioclavicular joint, starting over the clavicle. The incision should be curvilinear to the lines of Langer toward the coracoid process.
2. Uncoil the Nitinol wire and pass the looped end into the back of the instrument handle. Push the wire forward, guiding it through the opening at the thumb grip, until the loop sits just inside the tip of the hook.
3. Expose the superior surface of the coracoid and pass the tip of the hooked end underneath and around it from medial to lateral, following the surface of the bone until the tip of the hooked end visibly protrudes from the soft tissue.
4. Advance the Nitinol wire either at the thumb grip or at the back of the handle until the loop protrudes from the hooked end of the instrument.
5. Pass the end of the implant through the loop in the Nitinol wire and secure it by retracting the wire until it is trapped against the end of the instrument. Hold the Nitinol wire at the thumb grip to prevent the tape from slipping.
6. Carefully withdraw the instrument back around the coracoid process and out on the medial side, pulling the implant around with it.
7. Release the implant from the loop in the Nitinol wire and remove the instrument.

## Disposal

No specific disposal requirements other than handling contaminated items as sharps clinical waste.

## Complaints

Any health care professional who has any complaints or experienced any dissatisfaction in the product quality, identity, durability, reliability, safety, usability, effectiveness, and/or performance, should notify the manufacturer and distributor immediately.

If the product ever malfunctions and may have caused or

contributed to the death or serious injury of a patient, the manufacturer and relevant local regulatory authority should be notified immediately by telephone, email or written correspondence.

When filing a complaint, provide the component(s) name and number, lot number(s), your name and contact details and the nature of the complaint.

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All rights reserved. © Neoligaments™ 2020.
Neoligaments and Xiros are trademarks of
Xiros™ Ltd
Xiros Limited, registered in England No. 1664824.