

# Infinity-Lock™ Button System

- EN** Instructions for Use
- ES** Instrucciones de uso
- DE** Hinweise zur Anwendung
- EL** Οδηγίες χρήσης
- FR** Consignes d'utilisation
- IT** Istruzioni per l'uso
- PT** Instruções de uso
- TR** Kullanım Talimatları
- NL** Gebruiksaanwijzing
- SV** Användningsinstruktioner
- PL** Instrukcja użytkowania
- MS** Arahan penggunaan
- FI** Käyttöohjeet
- TH** คำแนะนำในการใช้
- KO** 사용 설명서


[1]


[2]


[3]


[4]


[5]


[6]


[7]

[8]

### Ordering Information (supplied sterile)

102-1089 Infinity-Lock Button System:

Infinity-Lock Tube-Tape, 7 mm x 240 mm

Infinity-Lock Button, 4 mm x 12 mm

Cannulated drill bit, plain shank to fit Jacobs chuck,  
4.0 mm diameter x 120 mm

Guidewire, diameter 2.0 mm x 150 mm

Optional Disposable Instrument:

202-1411 CC-Hook, with a curved end, Left

202-1413 CC-Hook, with a curved end, Right

### Description

The Infinity-Lock Button System comprises a permanent implantable 240 mm Tube-Tape and titanium alloy Button together with a disposable cannulated drill bit and guidewire. A coracoid passer, such as the Xiros CC-Hook, is also required.

### Material Specifications

The Tube-Tape is made from polyethylene terephthalate (polyester) and the Button from implant grade titanium alloy (Ti-6Al-4V) to ISO 5832-3. The cannulated drill and guidewire are made from stainless steel 17-4PH (ASTM F899) and 316 (ASTM A262) respectively. These materials meet applicable specifications established by ISO 10993-1 Biological evaluation of medical devices and are appropriate for their intended use.

### Intended Use

The Infinity-Lock Button System is intended to provide fixation during the healing process following a syndesmotic trauma, such as fixation of acromioclavicular separations due to coracoclavicular ligament disruption.

### Indications

The Infinity-Lock Button System is indicated for patients with acromioclavicular separations resulting from disruption to the coracoclavicular ligaments.

### Contraindications

- Known hypersensitivity to implant materials. If the patient is suspected of having any foreign body sensitivity, appropriate tests should be made prior to implantation.

- Infections or any structural or pathological condition of the bone or soft tissue that would be expected to impair healing or secure fixation.
- Patients unable or unwilling to restrict activities to prescribed levels or follow a rehabilitation programme during the healing period.
- Skeletally immature patients are not suited as the Infinity-Lock will not elongate with growth. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.

### Warnings

- The Infinity-Lock Button System is provided sterile and remains so unless the packaging is damaged or opened. Store in standard conditions. DO NOT RSTERILIZE.
- The Infinity-Lock Button System is for single use only as it would be damaged by reprocessing and the instruments are subject to wear. Do not use after the expiration date. Discard any open, unused product.
- The surgeon must be thoroughly familiar with these instructions and the surgical technique recommended overleaf before using the Infinity-Lock Button System.
- The general principles of patient selection and sound surgical judgement apply to the reconstruction procedure.
- The Infinity-Lock Button System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artefact in the MR environment. The safety of the Infinity-Lock Button System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.
- Revision cases where the recommended position of the bone tunnel has been compromised. In such cases there is a risk that the more oblique direction of the repair could cause the device to fail.

### Precautions

- Inspect the device, packaging and labelling prior to use and do not use if damaged. Contact Xiros if the package has been opened or altered.
- Care should be taken to avoid nerves and other anatomical hazards during surgery, as well as possible infection.
- Take care to round the upper edges of the bone tunnel to prevent abrasion of the Tube-Tape.
- Ensure there is sufficient bone around the tunnel in the clavicle to resist expected forces, taking account of the quality of the bone.
- Check that the repair is physiological and does not restrict range of motion.

- When trimming the Tube-Tape to length, cut straight across at right angles to minimize the generation of loose fibres and carefully remove any created from the wound.
- Avoid damage when handling the Tube-Tape. Avoid crushing or crimping when using surgical instruments such as forceps or needle holders.
- It is recommended that patients are placed on prophylactic antibiotics prior to surgery, to minimize the risk of latent infections developing at the implant site.
- Use aseptic technique throughout the procedure.
- Patients should be warned not to exceed appropriate activity levels or to overload the repair before complete healing has occurred.
- When drilling the bone tunnel ensure the tape is positioned away from the underside of the clavicle to prevent accidental damage from the drill.

### Potential Adverse Effects

- As with any surgical procedure there is a risk of infection.
- Allergic and inflammatory tissue reactions.

### Surgical Technique

1. With the patient in the beach chair position use a vertical 5 cm skin incision starting at the level of the clavicle and slightly medial to the tip of the coracoid. Incise the fascia and the deltoid vertically and then divide the periosteum over the posterior clavicle laterally as far as the AC joint.
2. Carry out a sub-periosteal dissection creating an "L" shaped flap then insert a stay suture into the apex of the flap to aid retraction.
3. Gently apply a self-retaining retractor to aid access to the coracoid process.
4. Insert a cholecystectomy-type forceps under the acromion and follow this with the coracoid passer, passed under the neck of the coracoid from medial to lateral. See the coracoid passer instructions for use for further details.
5. Capture the green lead suture of the Tube-Tape with the coracoid passer (Figure 1) and pull it medially under the coracoid so that the loop of the Tube-Tape is accessible (Figure 2). Remove the lead suture when satisfied.
6. Pass both limbs of the Tube-Tape through the loop, lassoing the coracoid, and pull vertically upwards.
7. Identify the tuberosity on the inferior surface of the clavicle to which the conoid ligament was attached before being avulsed.
8. Drill the 2 mm guidewire perpendicularly through the middle of the clavicle at the level of the conoid tuberosity.
9. Over-drill with the 4 mm cannulated drill bit

to create the final bone tunnel in the clavicle.

Chamfer the upper edges of the tunnel to prevent abrasion of the Tube-Tape. Take care to leave sufficient bone around the tunnel to resist expected forces.

10. Use a suture and needle to pass the limbs of the Tube-Tape through the bone tunnel one at a time. Alternatively this may be performed with the nitinol wire from the CC-Hook, if being used. (Figure 3).
11. Pass the limbs of the Tube-Tape through the central holes of the Infinity-Lock Button, using the suture and needle if necessary (Figure 4).
12. The acromial end of the clavicle is reduced to the normal anatomical position by pushing downward on it while the arm is simultaneously pulled upward.
13. While maintaining this reduction, apply appropriate tension to the Tube-Tape (Figure 4). Take care to avoid damaging the Tube-Tape. If desired an optional half-knot may be tied over the top of the Button (Figure 5).
14. Check the repair is physiological and does not affect range of motion. When satisfied pass the ends of the Tube-Tape either side of the clavicle and secure inferiorly with a surgeons knot (Figures 5, 6 & 7). If there is insufficient space below the clavicle the knot may be tied anteriorly. Ensure the knot is securely locked.
15. Cut any excess Tube-Tape with scissors at right angles to its length to minimize the generation of loose fibres. Take care to remove any loose fibres that are created. Stitch the cut ends of Tube-Tape back on itself (Figure 8).
16. Repair the soft tissues by re-attaching the "L" shaped flap while tensioning the superior acromioclavicular ligament during the repair. Ensure the cut ends of the Tube-Tape are well buried in tissue.

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All rights reserved. © Xiros 2018. Worldwide
Xiros are trademarks of Xiros.