Rota-Lok™ System Implant Set

EN Instructions for Use

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


[4]

Ordering Information (supplied sterile)

102-1064 Rota-Lok System Implant Set, includes:

Rota-Lok, 10 mm x 500 mm

Packaged with the following disposable:

Drill Bit, plain shank to fit Jacobs Chuck, 2.5 mm diameter

Description

The Rota-Lok is a permanent implantable device for reconstruction of chronic, massive, full thickness tears of the rotator cuff without the need to harvest any tissue. The implant set consists of a 10 mm x 500 mm Rota-Lok device and the required 2.5 mm drill bit. The Rota-Lok is looped through the torn rotator cuff then threaded through two small bone tunnels in the humeral head and secured by knotting. The open weave section of the Rota-Lok partially covers the superior aspect of the humeral head and bridges the gap to the torn rotator cuff, thereby restoring the humeral head to its normal position to maximize function of the deltoid.

Material Specifications

The Rota-Lok is made from polyethylene terephthalate (polyester) and the drill from instrument grade stainless steel to ASTM F899. These materials meet applicable specifications established by ISO 10993-1 Biological evaluation of medical devices and are appropriate for their intended use.

Intended Use

The Rota-Lok System is intended to be used in the reconstruction of chronic, massive, full thickness tears of the rotator cuff.

Indications

The Rota-Lok System is indicated for patients with retracted rotator cuff tears that cannot be mobilized back to the bony attachment site and where the tissue has undergone degeneration.

Contraindications

- Known hypersensitivity to materials. If the patient is suspected of having any foreign body sensitivity, appropriate tests should be made prior to implantation.

- Patients with cuff tear arthropathy, primary glenohumeral osteoarthritis with stiffness, or previous history of dislocation with recurrent instability.
- Infections or any structural or pathological condition of the bone or soft tissue that would be expected to impair healing or secure fixation.
- Patients unable or unwilling to restrict activities to prescribed levels or follow a rehabilitation programme during the healing period.
- Skeletally immature patients. The Rota-Lok will not elongate with growth and must not bridge or disrupt the growth plate.

Warnings

- The devices are provided sterile and remain so unless the packaging is damaged or opened. DO NOT RESTERILIZE the Rota-Lok.
- The devices are for single use only as they would be damaged by reprocessing. Do not use after the expiration date. Discard any open, unused product.
- The surgeon must be thoroughly familiar with these instructions and the recommended surgical procedure before using the devices.
- The general principles of patient selection and sound surgical judgement apply to the reconstruction procedure.

Precautions

- Inspect the devices, packaging and labelling prior to use and do not use if damaged. Contact the Neoligaments Sales Department if the packaging has been opened or altered. Store in standard conditions.
- Care should be taken to avoid the axillary nerve and other anatomical hazards during surgery.
- Ensure there is sufficient bone around the tunnels to resist expected forces, taking account of bone quality.
- Round bone tunnels where possible to prevent abrasion of the Rota-Lok.
- Check that the repair is physiological and the knot is securely locked and well buried in tissue.
- When trimming the Rota-Lok, cut straight across at right angles to minimize the generation of loose fibres and carefully remove any created.
- Avoid damage when handling the Rota-Lok. Avoid crushing or crimping when using surgical instruments such as forceps or needle holders.
- Aseptic technique should be followed throughout and wherever possible patients should be placed on prophylactic antibiotics prior to surgery to minimize the risk of latent infections developing at the implant site.
- Patients should be warned not to exceed appropriate activity levels before complete healing has occurred.

Potential Adverse Effects

- Wound dehiscence, infected wounds, minimal acute inflammatory tissue reaction and transitory local irritation.
- The Rota-Lok could possibly break through the bone if bone tunnels are incorrectly positioned, the bone is of poor quality, or prescribed activity levels are exceeded.
- The Rota-Lok integrates well with the patient's tissue and consequently may not be suitable where eventual removal of the Rota-Lok is anticipated.

Instructions for Use

The following provides instructions on how to use the Rota-Lok System during a standard approach.

1. The procedure is performed in the beach chair position under general anaesthesia, and/or interscalene block.
2. Make a lateral mini-open incision from the anterolateral margin of the acromion, extending approximately 5 cm inferiorly. Alternatively extend a standard lateral portal superiorly and inferiorly if evaluating arthroscopically first.
3. Dissect the skin and subcutaneous tissue, secure a haemostasis, split the deltoid along its fibres and insert a self-retaining retractor. Take care to avoid splitting the deltoid too far anteriorly, to avoid damage to the axillary nerve.
4. Apply a second retractor to separate the head from the acromion to visualize the cuff. Grasp the cuff with Kocher artery forceps and mobilize by blunt dissection.
5. Make a stab wound through a modified Neviaser (supraspinatus) portal, 1 cm medial to the apex of the posterior border of the clavicle and the medial border of the acromion to minimize the risk of disrupting nerves and arteries in this region.
6. Tension the cuff and pass a pair of Fraser Kelly artery forceps through the healthy posterior portion, sufficiently medial to the tear to prevent the Rota-Lok from pulling through the soft tissue. Avoid excessive force when pushing the forceps through the cuff to prevent damage to the underlying cartilage surface on the humeral head.
7. Grasp the corded end of the Rota-Lok with the Fraser Kelly artery forceps (Figure 1). Pull it through the cuff and out of the modified Neviaser portal until approximately 5 cm of the central open weave section remains protruding from the inferior surface of the cuff.
8. Introduce a curved hook over the upper surface of the cuff and draw the superior end of the Rota-Lok out of the lateral wound (Figure 2).
9. Repeat the steps above to pass the lower end of the Rota-Lok through the inferior surface of the cuff to its superior surface and then out of the lateral wound.
10. Lay both ends of the Rota-Lok flat on the superior surface of the cuff, ensuring that each has a similar length of open weave mesh protruding from the cuff.
11. Tension both ends of the Rota-Lok to seat the mesh tightly against the cuff. Pull the ends towards the original attachment site to mobilize the cuff.
12. Identify the location for two bone tunnels within the footprint of the greater tuberosity. For an anterosuperior tear one tunnel will typically be located in the anterior region and one slightly anterior to the middle. For a posterosuperior tear one tunnel will typically be located in the posterior region and one slightly posterior to the middle.
13. Use the supplied 2.5 mm drill bit to create two curved tunnels in the humeral head as described below. Leave sufficient bone around the tunnels to resist expected forceps, taking account of bone quality (Figure 3).
14. When drilling the anterior tunnel externally rotate the arm and drill a hole slightly obliquely. Mark the tunnel entrance with diathermy and insert a bone awl.
15. Drill an oblique hole at the outer side of the humerus to meet the tip of the bone awl and mark the tunnel entrance with diathermy.
16. Repeat these steps with the arm internally rotated to create a posterior tunnel parallel to the first.
17. Use a bone awl to ream the intersections of the two parallel tunnels to ensure they are smooth and free of obstruction. Where possible round the tunnel edges to prevent abrasion of the Rota-Lok.
18. Whilst maintaining tension on the rotator cuff use either a USP #2 suture with a curved needle or a suture snare to pass each corded end of the Rota-Lok through the tunnel nearest to it.
19. Hold the arm in 20-30° abduction and tension the cords to ensure the length of the reconstruction is physiological and covers at least the superior convex surface of the humeral head. Avoid over- or under-tensioning to ensure adequate joint function can be achieved.
20. Whilst maintaining tension check stability and range of motion before tying the corded ends of the Rota-Lok on the outside of the humerus using a reef knot with an extra throw (Figure 4).
21. Recheck shoulder movement and ensure the knot is locked before using scissors to trim the cords to length. Cut each cord at right angles to its length to minimize the generation of loose fibres, taking care to remove any that are created.
22. Take care to ensure the knot remains buried

in tissue. If necessary the cut ends may be restrained by stitching them back to the cord.

23. Suture the remnant of the cuff to the flattened Rota-Lok to close the defect as much as possible.

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All rights reserved. © Neoligaments™ 2017.
trademarks of Xiros.
Xiros Limited, Registered in England No. 1664824.