AchilloCord^{PLUS}™ System Implant Set

EN Instructions for Use

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

AchilloCord^{PLUS} System Implant Set (Product Code 102-1142) includes:

AchilloCord^{PLUS}, 5 mm x 800 mm

Packaged with the following disposables:
Rigid Probe with eye, stainless steel, 20 cm
Drill Bit plain shank to fit Jacobs Chuck,
3.2 mm diameter

Patient Information

The following information is provided for use by clinicians, however as the learned intermediary between the company and the patient, the clinician must convey the aspects they consider relevant to the individual patient. The patient must be informed of the potential adverse effects (risks/complications) contained in this insert (see **POTENTIAL ADVERSE EFFECTS**).

Description

The AchilloCord^{PLUS} is a permanent implantable device for reconstruction of an acute Achilles tendon rupture. It is 5 mm in diameter by 800 mm in length and has a densely woven, flexible tubular structure that spreads out as it is loaded to prevent cutting through the tissue. It is passed through a bone tunnel in the os calcis and the tendon stumps before being knotted in a subcutaneous Bunnell-type suture. The AchilloCord^{PLUS} System is packaged with the required 3.2 mm drill bit and 20 cm rigid probe. The AchilloCord^{PLUS} is supplied for single use only and is intended for implantation by surgeons qualified in acute Achilles tendon surgical repair.

Material Specifications

The AchilloCord^{PLUS} is made from 100% polyethylene terephthalate (polyester).

Intended Use

The AchilloCord^{PLUS} is intended to be used in the reconstruction of acute Achilles tendon ruptures.

Indications for Use

The AchilloCord^{PLUS} is indicated for patients with acute ruptures of the Achilles tendon but is particularly suited to patients where an extended period of postoperative immobilization is undesirable.

Contraindications

The AchilloCord^{PLUS} should not be implanted in patients:

- With known hypersensitivity to implant materials. If the patient is suspected of having any foreign body sensitivity, appropriate tests should be made prior to implantation.
- With infections or any structural or pathological condition of the bone or soft tissue that would be expected to impair healing or prevent secure fixation.
- With chronic ruptures of the Achilles tendon.
- That are skeletally immature as the AchilloCord^{PLUS} will not elongate with growth.
- Unable or unwilling to restrict activities to prescribed levels or to follow a rehabilitation program during the healing period.

Warnings

Pre-operative

- The device is for single use only as reprocessing would damage it or could lead to infection.
- Do not use after the expiration date.
- Users must be thoroughly familiar with these instructions and the recommended surgical procedure before using the device.
- The general principles of patient selection and sound surgical judgement apply to the reconstruction procedure.

MRI Safety information

The AchilloCord^{PLUS} is MR safe (i.e. an item that poses no known hazards) in all MR environments.

Precautions

Packaging

- The device is provided sterile (by gamma irradiation) and remains so unless the packaging is damaged or opened.
- Do not re-sterilize.
- Inspect the device, packaging and labelling prior to use and do not use if damaged. Contact the Neoligaments Sales Department if the package has been opened or altered. Damaged packages or products should be returned to Neoligaments.

Handling and storage

- Store in standard conditions.
- Discard any open, unused product.
- All devices should be inspected prior to use to ensure that there is no deterioration such as discoloration or damage to the fibres.
- Avoid damage when handling the AchilloCord^{PLUS}. Avoid crushing or crimping when using surgical instruments such as forceps or needle holders.

Pre-operative

The safety and effectiveness of this device has not been established in the following patient populations:

- Pregnant and breastfeeding women.

- It is recommended that wherever possible patients are placed on prophylactic antibiotics prior to surgery, to minimize the risk of infection developing at the implant site.

Intra-operative

- When passing the AchilloCord^{PLUS} back through a stab incision, pull it to confirm there is no puckering of the skin that indicates trapped subcutaneous tissue. Loosen the tissue with artery forceps if required.
- Care should be taken to avoid the sural nerve during surgery, as well as possible infection.
- Ensure there is sufficient bone around the tunnel to resist expected forces, taking account of the quality of the bone.
- Check that the repair is physiological, and that the knot is securely locked and well buried.
- When trimming the AchilloCord^{PLUS}, cut straight across at right angles to minimize the generation of loose fibres and carefully remove any created from the wound.

Post-operative

- The patient should be warned not to exceed appropriate activity levels before complete healing has occurred.

Potential Adverse Effects

Below is a list of the potential adverse effects (e.g., complications) associated with the use of the device. Additional surgery may be required to correct some of these events.

As with any procedure of this type, there is a risk that surgery may not be effective in treatment or may cause worsening symptoms.

- Pertinent risks associated with any surgical procedure include: Wound dehiscence, superficially infected wounds and deep infection.
- Pertinent risks associated with acute Achilles tendon surgical repair include: Sural nerve damage, Achilles tendon re-rupture, deep vein thrombosis.
- Risks associated with synthetic ligaments, including the AchilloCord^{PLUS}, include: Implant rupture, scar tethering, minimal acute inflammatory tissue reaction, allergic reaction and transitory local irritation.

Surgical Technique

The following provides instructions on how to use the AchilloCord^{PLUS} during a standard approach.

- Make a 2.5-3 cm posteromedial incision at the site of the rupture, parallel to the medial border of the Achilles tendon. Dissect the underlying soft tissues to expose the site of the injury. Take care to avoid damaging the sural nerve.
- Make a stab incision on the lateral side of the os calcis and use the supplied drill to create a transverse tunnel along the long axis of the tendon

approximately 1-2 cm below its insertion. Take care not to compromise the strength of the tunnel by drilling too close to the bone surface. As the drill emerges from the medial side make a stab incision to allow it to exit.

- Thread the AchilloCord^{PLUS} through the eyelet of the supplied probe and pass it through the bone tunnel from the lateral side until approximately 25 cm emerges.
- Use the probe to make a subcutaneous path from the medial stab incision, running as close as possible to the bone surface then through the distal tendon stump, centrally in the sagittal plane, exiting at the centre of the rupture (Figure 1). Pull the AchilloCord^{PLUS} through so that it disappears under the skin on the medial wall of the os calcis. Use artery forceps to release any trapped subcutaneous tissue if necessary.
- Remove the probe and rethread onto the lateral end of the AchilloCord^{PLUS}. Make another proximal subcutaneous path through the distal tendon stump from the lateral side in the same manner as above. Tension both ends once they have exited the distal tendon stump to remove any slack.
- If possible approximate the ends of the Achilles tendon and pass the probe through the centre of the proximal stump, using a stab incision if necessary, so that it emerges about 2 cm proximal to the site of injury. Create a stab incision on the medial side of the leg where the probe tents the skin and pull the probe and AchilloCord^{PLUS} through, while maintaining tension on the distal tendon stump.
- Pass the probe back through the same medial stab incision to the lateral side, exiting the skin through a stab incision about 4 cm proximal to the site of injury. Avoid damaging the sural nerve or blood supply in the muscle or tendon.
- Pass the probe horizontally back through the same lateral stab incision to exit a stab incision at the medial side. Avoid trapping the sural nerve as the AchilloCord^{PLUS} is looped back through the Achilles tendon.
- Pass the probe back through the same medial stab incision, exiting at a lateral stab incision about 2 cm proximal to the site of injury.
- Make the final pass with the probe distally through the same lateral stab incision to exit at the proximal tendon stump. Exit at the deep portion of the stump to help place the knot away from the skin.
- Pull the ends of the AchilloCord^{PLUS} to ensure there is no slackness within the tendon and confirm there is no trapped subcutaneous tissue by checking for puckering of the skin at the stab incisions. Release the tissue with artery forceps if necessary.
- Place the ankle in 5-10° plantar flexion and tension the proximal end of the AchilloCord^{PLUS} to

remove any slack and approximate the two ends of the Achilles tendon. Pull the distal end of the AchilloCord^{PLUS} to remove any bunching of the cord at the site of injury.

13. Tie the ends of the AchilloCord^{PLUS} at the site of injury with either a surgeon's knot or a reef knot, with or without an extra throw as required. Ensure the knot is locked (Figure 2).
14. Avoid over-/under-tightening or shortening of the tendon to ensure the length of the reconstruction is physiological and the AchilloCord^{PLUS} is appropriately taut.
15. Evaluate the position of the tendon when the knee is flexed to 90° and in full extension. Assess the repair with Simmonds' test by squeezing the calf muscle to demonstrate plantar flexion of the foot.
16. If satisfied cut any excess AchilloCord^{PLUS} with scissors at right angles to its length to minimize the generation of loose fibres. Remove any loose fibres that are created and suture the cut ends back to the cord.
17. Use a strong non-absorbable suture to stitch the ends of the tendon together over the top of the knot so that it is well buried in the tissue.

Disposal

No specific disposal requirements other than handling contaminated items as clinical waste.

Device Retrieval Methods

Should it be necessary to remove an AchilloCord^{PLUS}, contact Neoligaments prior to the scheduled surgery for product/tissue retrieval information.

Complaints

Any health care professional who has any complaints or experienced any dissatisfaction in the product quality, identity, durability, reliability, safety, usability, effectiveness, and/or performance, should notify the manufacturer and distributor immediately.

If the product ever malfunctions and may have caused or contributed to the death or serious injury of a patient, the manufacturer and relevant local regulatory authority should be notified immediately by telephone, email or written correspondence.

When filing a complaint, provide the component(s) name and number, lot number(s), your name and contact details and the nature of the complaint.

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All rights reserved. © Neoligaments™ 2020.
Neoligaments and Xiros are trademarks of
Xiros Limited, Registered in England No. 1664824.