AchilloTape™ System

EN Instructions for Use

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
Ordering Information (supplied sterile)

Small

102-1076 AchilloTape Procedure Pack (10 mm x 800 mm), includes:

Poly-Tape, 10 mm x 800 mm

Packaged with the following disposables:

Rigid Probe with eye – stainless steel, 20 cm
Drill Bit, plain shank to fit Jacobs Chuck, 3.2 mm diameter

Large

102-1077 AchilloTape Procedure Pack (20 mm x 800 mm), includes:

Poly-Tape, 20 mm x 800 mm

Packaged with the following disposables:

Rigid Probe with eye – stainless steel, 20 cm
Drill Bit, plain shank to fit Jacobs Chuck, 3.6 mm diameter

Description

The AchilloTape is a permanent implantable device for repair of Achilles tendon ruptures. The procedure pack contains an 800 mm long Poly-Tape with an open weave polyester mesh, a 20 cm probe and a drill bit. The procedure pack is available in small and large sizes with the Poly-Tape 10 mm or 20 mm wide and the drill bits 3.2 mm or 3.6 mm in diameter respectively. The repair is performed by using the probe to pass the Poly-Tape through a bone tunnel drilled in the os calcis then the ruptured tendon stumps before being knotted in a Bunnell type suture.

Material Specifications

The Poly-Tape is made from polyethylene terephthalate (polyester) and the probe and drills from instrument grade stainless steel (ASTM F899). These materials meet applicable specifications established by

ISO10993-1 Biological evaluation of medical devices and are found to be appropriate for their intended use.

Intended Use

The AchilloTape is intended for repair of a ruptured Achilles tendon.

Indications

The AchilloTape is indicated for patients with ruptures of the Achilles tendon but is particularly suited to the following cases:

- Delayed diagnosis of chronic ruptures with tendon atrophy.
- Previous failed surgery or conservative treatment.
- Undesirable prolonged immobilization for acute ruptures.
- Patients with a history of prolonged steroid treatment or the elderly where wound healing is often a problem.

Contraindications

- Known hypersensitivity to implant or instrument materials. If the patient is suspected of having any foreign body sensitivity, appropriate tests should be made prior to implantation.
- Infections or any structural or pathological condition of the bone or soft tissue that would be expected to impair healing or secure fixation.
- Patients unable or unwilling to restrict activities to prescribed levels or follow a rehabilitation programme during the healing period.
- Poly-Tapes integrate well with soft tissue so may not be suitable where eventual removal of the device is anticipated.
- The Poly-Tape is not suitable for skeletally immature patients as it will not elongate with growth.

Warnings

- The devices are provided sterile and remain so unless the packaging is damaged or opened. Store in standard conditions. DO NOT RESTERILIZE the Poly-Tape
- The devices are for single use only as they would be damaged by reprocessing. Do not use after the expiration date. Discard any open, unused product.
- The surgeon must be thoroughly familiar with these instructions and the recommended surgical procedure before using the devices.
- The general principles of patient selection and sound surgical judgement apply to the reconstruction procedure.
- Select the most suitable size of AchilloTape for the surgery to be undertaken, based on the patient's anatomy, build and lifestyle, size of rupture, or extent of poor quality tissue.

Precautions

- Inspect the device, packaging and labelling prior to use and do not use if damaged. Contact the Neoligaments sales department if the package has been opened or altered.
- When passing the Poly-Tape back through a stab incision pull it to confirm there is no puckering of the skin that indicates trapped subcutaneous tissue. Loosen the tissue with artery forceps if required.
- Care should be taken to avoid the sural nerve and other anatomical hazards during surgery, as well as possible infection.
- Ensure there is sufficient bone around the tunnel in the tibial tuberosity to resist expected forces, taking account of the quality of the bone. Round any bone tunnel edges where possible.
- Check that the repair is physiological and the knot is securely locked and well buried in tissue.
- When trimming the Poly-Tape cut straight across at right angles to minimise the generation of loose fibres and carefully remove any created from the wound.
- Avoid damage when handling the Poly-Tape. Avoid crushing or crimping when using surgical instruments such as forceps or needle holders.
- It is recommended that wherever possible patients are placed on prophylactic antibiotics prior to surgery, to minimize the risk of latent infections developing at the implant site.
- Patients should be warned not to exceed appropriate activity levels or to overload the repair before complete healing has occurred.
- The patient should be warned to protect the wound from the environment to prevent infection and poor wound healing.
- It is advisable to use a plaster back slab to protect the repair.

Potential Adverse Effects

- Wound dehiscence, infected wounds, minimal acute inflammatory tissue reaction, and transitory local irritation.
- The Poly-Tape could possibly break through the bone if bone tunnels are incorrectly positioned, the bone is of poor quality, or prescribed activity levels are exceeded.

Surgical Technique

Chronic Repair

1. Make a 6-9 cm posteromedial incision at the site of the rupture, extending proximally more than distally (Figure 1). Dissect the underlying soft tissues to expose the site of the injury then mobilize and approximate the tendon stumps.

Take care to avoid damaging the sural nerve.

2. Make a stab incision on the lateral side of the os calcis and using the supplied drill bit, drill a transverse tunnel along the long axis of the tendon approximately 1-2 cm below its insertion, taking care to leave sufficient bone to resist expected forces. As the drill emerges from the medial side make a stab incision to allow it to exit (Figure 1).
3. Select the most suitable size of AchilloTape for the surgery to be undertaken based on the patient's anatomy, build and lifestyle, size of rupture, or extent of poor quality tissue.
4. Thread the Poly-Tape through the eyelet of the supplied probe. Pass it through the bone tunnel from the lateral side until approximately 25 cm emerges.
5. Use the probe to make a subcutaneous path from the medial stab incision, running as close as possible to the bone surface then through the distal tendon stump, centrally in the sagittal plane, exiting at the centre of the rupture.
6. Pull the Poly-Tape through so that it disappears under the skin on the medial wall of the os calcis. Use artery forceps to release any trapped tissue if necessary.
7. Remove the probe and rethread onto the lateral end of the Poly-Tape. Make a subcutaneous path through the distal tendon stump from the lateral side in the previous manner. Cross the ends of the Poly-Tape at the exit of the tendon stump and tension to remove any slack (Figure 1).
8. Pull the proximal end of the tendon distally and keep under tension. Ideally the gap should close but may remain if necessary to achieve correct biomechanics of the tendon.
9. Use the probe to pass the long end of the Poly-Tape through the centre of the proximal stump so that it emerges medially about 2 cm proximal to the site of injury. Use a stab incision in the centre of the tendon if necessary.
10. Follow the steps below to create a Bunnell-type repair in the tendon while minimising the distance between entry and exit points to reduce the amount of Poly-Tape on the tendon surface.
11. Use the probe to pass the Poly-Tape through the tendon at an angle to exit on the lateral side about 4 cm proximal to the site of injury. Avoid damaging the sural nerve or blood supply in the muscle or tendon.
12. Pass the probe horizontally back through the tendon close to the exit at the medial side. Avoid trapping the sural nerve as the Poly-Tape is looped through the Achilles tendon.
13. Pass the probe back through the tendon at an angle, exiting on the lateral side about 2 cm proximal to the site of injury.

-
14. Make the final pass with the probe distally to exit at the deep portion of the proximal tendon stump to help place the knot away from the skin.
 15. Pull the ends of the Poly-Tape to ensure there is no slackness within the tendon and confirm there is no trapped subcutaneous tissue by checking for puckering of the skin around the stab incisions. Release the tissue with artery forceps if necessary.
 16. Place the ankle in 5-10° plantar flexion and tension the proximal end of the Poly-Tape to remove any slackness and approximate the two ends of the Achilles tendon. Pull the distal end of the Poly-Tape to remove any bunching of the tape at the site of rupture.
 17. Tie the ends of the Poly-Tape at the site of injury with either a surgeon's knot or a reef knot, with or without an extra throw as required. Ensure the knot is locked.
 18. Avoid over- or under-tightening or shortening of the tendon to ensure the length of the reconstruction is physiological and the Poly-Tape is appropriately taut.
 19. Evaluate the position of the tendon when the knee is flexed to 90° and in full extension. Assess the repair with the Simmonds' test by squeezing the calf muscle to demonstrate plantar flexion of the foot.
 20. If satisfied cut any excess Poly-Tape with scissors at right angles to its length to minimize the generation of loose fibres. Remove any loose fibres that are created and suture the cut ends back to the Poly-Tape.
 21. If the ends of the ruptured tendon can be approximated, bury the knot deep to the tendon (Figure 2) and use a strong non-absorbable suture to stitch the ends of the tendon together over the top of the knot.
 22. If the ends of the tendon cannot be approximated and a gap remains, mobilize and split the scar tissue then bring it over and suture it to minimize the amount of Poly-Tape and knot exposed.

Acute Repair

When performing an acute repair, follow the description for the chronic repair with the following changes:

In step 1 use a small posteromedial incision 2.5-3 cm long.

In step 7 it is not essential to cross the strands of the Poly-Tape in the distal tendon stump.

In steps 9-13 pass the Poly-Tape through stab incisions where the probe tents the skin to create the Bunnell-type repair (Figure 3).

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All rights reserved. © Neoligaments™ 2017.
trademarks of Xiros.
Xiros Limited, Registered in England No. 1664824.