

Ortho-Tape

High Strength Narrow Sterile Polyester Tape

Benefits

Made from woven high tenacity polyester, providing high strength and high stiffness

Heat sealed ends prevent fraying and facilitate easier passage through bone tunnels and soft tissue

Narrower sized 3mm tape provides higher strength than competitor 5mm wide braided devices

Indications

Ortho-Tape is intended for tissue approximation, including use in reconstructing damaged or torn ligaments, tendons and other soft tissues according to the surgeon's own preferred technique, and at his/her discretion.




Chart comparing ultimate tensile strength (N) and overall stiffness (N/mm) of the 3mm Ortho-Tape and two other brands.

Ordering Information

Multipack of 5 Ortho-Tapes (supplied sterile)

Ref	Weave	Width	Length	Strength (N)	
102-1028	Dense	3mm	400mm	690	
102-1027	Dense	3mm	600mm	690	

1 Ortho-Tape (supplied sterile with an 8mm Fastlok™)

Ref	Weave	Width	Length	
102-1331	Dense	3mm	400mm	

2 Ortho-Tapes (supplied sterile with a 6mm Fastlok™)

Ref	Weave	Width	Length	
102-1354	Dense	3mm	400mm	

For details of other Neoligaments polyester tapes, contact our Customer Services Department or visit our website at www.neoligaments.com

All rights reserved. © Neoligaments™ 2009. Worldwide patents and patents pending. Neoligaments, Xiros and Fastlok are trademarks of Xiros.



Developed and manufactured by **Neoligaments™**
 A division of Xiros™
 Springfield House Whitehouse Lane Leeds LS19 7UE
 Tel. +44 (0) 113 238 7202 Fax. +44 (0) 113 238 7201
enquiries@neoligaments.com
www.neoligaments.com